

Datasheet

Vixarelimab Biosimilar - Anti-OSMR mAb - Research Grade PX-TA1737

➤ Overview

Antibody name:	Vixarelimab Biosimilar - Anti-OSMR mAb - Research Grade
Isotype:	IgG4-IgG1,Lambda
Host species:	Homo sapiens
Production Host:	XtenCHO
Source:	CAS 2243320-83-2
Aliases/Synonyms:	30C6 (N73D) IGG4P.AGLY/IGG1, MMUNOGLOBULIN G (234-PROLINE,303-GLUTAMINE), ANTI-(HUMAN ONCOSTATIN M RECEPTOR SUBUNIT .BETA.) (HUMAN MONOCLONAL KPL-716 .GAMMA.-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL KPL-716 .LAMBDA.-CHAIN, DIMER, BIIB069,BIIB-069

➤ Properties

Clonality:	Monoclonal Antibody
Final buffer:	PBS buffer PH7.5
Purity:	>95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.