

Datasheet

Vixtimotamab Biosimilar - Anti-CD3E;CD33 mAb - Research Grade PX-TA1738

➤ Overview

Antibody name:	Vixtimotamab Biosimilar - Anti-CD3E;CD33 mAb - Research Grade
Isotype:	Bispecific Homodimer (VK-VH-VL'-VH', Tandem diabody);Kappa;Lambda
Host species:	Chimeric
Production Host:	XtenCHO
Source:	CAS 2243775-32-6
Aliases/Synonyms:	AMV564, MABFRAG HUMAN (FAB) ANTI P20138 (CD33_HUMAN)MABFRAG HUMAN (FAB) ANTI P07766 (CD3_HUMAN) (AMV564), CNTO 3953 TANDAB,T-652,CNTO-3953, AMV-564

➤ Properties

Clonality:	Monoclonal Antibody
Final buffer:	PBS buffer PH7.5
Purity:	>95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.