

## Datasheet

### Bavunalimab Biosimilar - Anti-LAG3;CTLA4 mAb - Research Grade PX-TA1813

#### ➤ Overview

**Antibody name:** Bavunalimab Biosimilar - Anti-LAG3;CTLA4 mAb - Research Grade  
**Isotype:** IgG1 Kappa;na Kappa  
**Host species:** Bispecific Mixed mAb and scFv  
**Production Host:** XtenCHO

#### ➤ Properties

**Clonality:** Monoclonal Antibody  
**Final buffer:** PBS buffer PH7.5  
**Purity:**

### Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

### Note:

For research use only. Not suitable for clinical or therapeutic use.