

## Datasheet

### Osemitamab Biosimilar - Anti-CLDN18 mAb - Research Grade PX-TA1872

#### ➤ Overview

|                         |  |
|-------------------------|--|
| <b>Antibody name:</b>   | Osemitamab Biosimilar - Anti-CLDN18 mAb - Research Grade |
| <b>Isotype:</b>         | IgG1 Kappa   |
| <b>Host species:</b>    | Homo Sapiens   |
| <b>Production Host:</b> | XtenCHO  |

#### ➤ Properties

|                      |                     |
|----------------------|---------------------|
| <b>Clonality:</b>    | Monoclonal Antibody |
| <b>Final buffer:</b> | PBS buffer PH7.5    |
| <b>Purity:</b>       |                     |

### Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

### Note:

For research use only. Not suitable for clinical or therapeutic use.