

Datasheet

e omega DT/Diphtheria toxin Antibody (DTD4) Biosimilar - Anti-Diphtheria toxin mAb - Research Grade PX-TA1903

➤ Overview

Antibody name: e omega DT/Diphtheria toxin Antibody (DTD4) Biosimilar - Anti-Diphtheria toxin mAb - Research Grade

Isotype: IgG2-kappa

Host species: Human

Production Host: XtenCHO

Aliases/Synonyms: 0

➤ Properties

Clonality: Monoclonal Antibody

Final buffer: PBS buffer PH7.5

Purity:

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.