

2'-O-Methyltransferase (NSP16)

PX-COV-P007

OVERVIEW

SIZE	100 ug
ORIGIN SPECIES	SARS-COV2
FRAGMENT	Full length
PROTEIN DELIVERED WITH TAG	N-terminal GST tag
MOLECULAR WEIGHT WITH TAG IF ANY	61kDa
DELIVERY CONDITION	Dry Ice


PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.Coli
PURITY	85%
PROTEIN ACCESSION	MN908947
FORM	Frozen
BUFFER	PBS, pH7.5, 0.02%NLS
STABILITY & STORAGE	Proteogenix did not perform stability or storage tests for your specific end products. Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

NSP16



Final sample QC. Coomassie blue staining

For research use only.