

3C-Like-protease

PX-COV-P003

OVERVIEW

SIZE	100ug
ORIGIN SPECIES	SARS-CoV2
FRAGMENT	Full length
PROTEIN DELIVERED WITH TAG	N-terminal Histidine tag
MOLECULAR WEIGHT WITH TAG IF ANY	34.95kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	80%
PROTEIN ACCESSION	MN908947
FORM	Frozen
BUFFER	PBS, pH7.5


STABILITY & STORAGE

ProteoGenix did not perform stability or storage tests for your specific end products. Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

3CL-PRO
3CLp
nsp5



Final sample QC. Coomassie blue staining

For research use only.