

## Recombinant NSP13 protein

PX-COV-P031

### OVERVIEW

---

<b>ORIGIN SPECIES</b>	SARS-CoV2
<b>FRAGMENT</b>	Full length
<b>PROTEIN DELIVERED WITH TAG</b>	N-terminal Histidine tag
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	69kDa
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	E.Coli
<b>PURITY</b>	>90%
<b>PROTEIN ACCESSION</b>	QHD43415.1
<b>FORM</b>	Frozen liquid
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

### More info

#### SYNONYMS

Non-structural protein 13

