

Recombinant NSP15 protein

PX-COV-P030

OVERVIEW

ORIGIN SPECIES	SARS-CoV2
FRAGMENT	Full length
PROTEIN DELIVERED WITH TAG	N-terminal Histidine tag
MOLECULAR WEIGHT WITH TAG IF ANY	41kDa
DELIVERY CONDITION	Dry Ice


PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.Coli
PURITY	>90%
PROTEIN ACCESSION	QHD43415.1
FORM	Frozen liquid
BUFFER	PBS, pH7.5
STABILITY & STORAGE	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

SARS-CoV-2 Non-structural protein-15;
Uridylate-specific
endoribonuclease



Final sample QC. Coomassie blue staining

For research use only.