

Recombinant NSP2 protein


PX-COV-P029

OVERVIEW

ORIGIN SPECIES	SARS-CoV2
FRAGMENT	Full length (Arg 1-Gly638)
PROTEIN DELIVERED WITH TAG	N-terminal Histidine tag
MOLECULAR WEIGHT WITH TAG IF ANY	72kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.Coli
PURITY	>90%
PROTEIN ACCESSION	YP_009725298.1
FORM	Frozen liquid
BUFFER	PBS, pH7.5
STABILITY & STORAGE	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles



More info

SYNONYMS

SARS-CoV-2 Non-structural protein - 2;
SARS-CoV-2 p65 homolog