

Recombinant RBD-Fc protein

PX-COV-P047

OVERVIEW

ORIGIN SPECIES	SARS-CoV2
FRAGMENT	RBD domain of Spike protein Arg319-Phe541
PROTEIN DELIVERED WITH TAG	C-terminal Fc tag
MOLECULAR WEIGHT WITH TAG IF ANY	53kDa
DELIVERY CONDITION	Dry Ice


PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	HEK293
PURITY	>90%
PROTEIN ACCESSION	YP_009724390.1
FORM	Frozen liquid
BUFFER	PBS, pH7.5
STABILITY & STORAGE	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

SARS-CoV-2 Receptor
Binding domain of Spike protein



Final sample QC. Coomassie blue staining

For research use only.