

RdRp catalytic domain of RNA-directed RNA polymerase

PX-COV-P006

OVERVIEW

ORIGIN SPECIES	SARS-CoV-2
FRAGMENT	Domain Asp4891-Val5212
PROTEIN DELIVERED WITH TAG	N-terminal GST tag
MOLECULAR WEIGHT WITH TAG IF ANY	64kDa
DELIVERY CONDITION	Dry Ice


PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.Coli
PURITY	>90%
PROTEIN ACCESSION	QHD43415.1
FORM	Frozen liquid
BUFFER	PBS, pH7.5, 0.02%NLS
STABILITY & STORAGE	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

SARS-CoV-2
Non-structural protein12



Final sample QC. Coomassie blue staining

For research use only.