

Recombinant human TMPRSS2 protein

PX-COV-P044

OVERVIEW

ORIGIN SPECIES	Human
FRAGMENT	Ser284-Gly492
PROTEIN DELIVERED WITH TAG	N-terminal Histidine tag
MOLECULAR WEIGHT WITH TAG IF ANY	24kDa
DELIVERY CONDITION	Dry Ice


PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.Coli
PURITY	>90%
PROTEIN ACCESSION	O15393
FORM	Frozen liquid
BUFFER	PBS, pH7.5, 0.02%NLS
STABILITY & STORAGE	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

Human transmembrane
Serine Protease 2



Final sample QC. Coomassie blue staining

For research use only.