

## Nucleoprotein

PX-COV-P010

### OVERVIEW

ORIGIN SPECIES	SARS-CoV2
FRAGMENT	Full length
PROTEIN DELIVERED WITH TAG	N-terminal Histidine tag + thrombin site
MOLECULAR WEIGHT WITH TAG IF ANY	47kDa
DELIVERY CONDITION	Dry Ice


### PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	85%
PROTEIN ACCESSION	MN908947
FORM	Liquid
BUFFER	PBS, pH7.5
STABILITY & STORAGE	Storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

### More

#### SYNONYMS

SARS-CoV-2  
Nucleocapsid protein



Final sample QC. Coomassie blue staining