

Recombinant S2 protein

PX-COV-P045

OVERVIEW

ORIGIN SPECIES	SARS-CoV2
FRAGMENT	Partial domain of Spike protein
PROTEIN DELIVERED WITH TAG	N-terminal Histidine tag
MOLECULAR WEIGHT WITH TAG IF ANY	60kDa
DELIVERY CONDITION	Dry Ice


PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	>90%
PROTEIN ACCESSION	QHD43416.1
FORM	Frozen liquid or lyophilized
BUFFER	PBS, pH7.5, 0.02% NLS
STABILITY & STORAGE	Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles

More info

SYNONYMS

Spike glycoprotein - S2 domain



Final sample QC. Coomassie blue staining

For research use only.