

## Anti-ECFP Mouse Monoclonal Antibody (8F6)

### OVERVIEW

---

PRODUCT REFERENCE	PTX1055
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1

### PRODUCT INFORMATION

---

EXPRESSION SYSTEM	Rabbit
BUFFER	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
STORAGE CONDITIONS	<p>Store at 4 °C for short term. Freeze and store at -80 °C for long term.</p> <p>Addition of 20 to 40% glycerol may improve cryoprotection.</p> <p>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.</p> <p>Aliquot the product according to your needs to avoid freeze/thaw cycles.</p>