



Anti- β -tubulin Mouse Monoclonal Antibody(Zebrafish Specific)

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PTX1092 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1 |

PRODUCT INFORMATION

| | |
|---------------------------|--|
| EXPRESSION SYSTEM | Mouse |
| BUFFER | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |

