



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti eGFP mouse monoclonal antibody

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PTX17860            |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | IgG1                |

### PRODUCT INFORMATION

---

|                           |  |
|---------------------------|--|
| <b>EXPRESSION SYSTEM</b>  | Mouse  |
| <b>BUFFER</b>             | PBS pH7.4, containing 0.02% NaN <sub>3</sub> and 50% glycerol  |
| <b>STORAGE CONDITIONS</b> | Store at 4 °C for short term. Freeze and store at -80 °C for long term.<br>Addition of 20 to 40% glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Aliquot the product according to your needs to avoid freeze/thaw cycles. |



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**