

Mouse IgA Isotype Control antibody (S107)

OVERVIEW

PRODUCT REFERENCE	PTX17906
CLONALITY	Monoclonal Antibody
ISOTYPE	IgA

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	PBS pH7.5
STORAGE CONDITIONS	<p>Store at 4 °C for short term. Freeze and store at -80 °C for long term.</p> <p>Addition of 20 to 40% glycerol may improve cryoprotection.</p> <p>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.</p> <p>Aliquot the product according to your needs to avoid freeze/thaw cycles.</p>