



## Mouse IgG2a LALA Isotype Control antibody (HyHEL-10)

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PTX17915            |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | IgG2A LALA          |

### PRODUCT INFORMATION

---

|                           |  |
|---------------------------|--|
| <b>EXPRESSION SYSTEM</b>  | Mouse  |
| <b>BUFFER</b>             | PBS pH7.5  |
| <b>STORAGE CONDITIONS</b> | Store at 4 °C for short term. Freeze and store at -80 °C for long term.<br>Addition of 20 to 40% glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Aliquot the product according to your needs to avoid freeze/thaw cycles. |

