



Mouse IgG2C Isotype Control antibody (HyHEL-10)

OVERVIEW

PRODUCT REFERENCE	PTX17924
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2C

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.

