



ANGPT1 Polyclonal Polyclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX18466
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	ANGPT1

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Human
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

