



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Goat Anti-Human IgG Polyclonal Antibody, FITC-conjugated

OVERVIEW

PRODUCT REFERENCE	PTX18862
CLONALITY	Polyclonal Antibody
TARGET	Anti-human antibodies

PRODUCT INFORMATION

HOST	Goat
REACTIVITY	Human
BUFFER & STORAGE	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.