



## Anti-Obinutuzumab Polyclonal Antibody

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX18871
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	Obinutuzumab

### PRODUCT INFORMATION

---

<b>HOST</b>	Rabbit
<b>REACTIVITY</b>	Drug
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4, 50% glycerol, 0.09% sodium azide.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

