



Anti-Human CD275/ICOSLG Nanobody (SAA1282)

OVERVIEW

PRODUCT REFERENCE	PTX18903
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD275/ICOSLG

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

