



# Proteogenix

## Proteogenix S.A.S

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## Proteogenix INC.

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |   |
|--------------------------|---|
| <b>PRODUCT NAME</b>      | Anti-SARS-CoV-2 S Protein VHH (SAA0974) |
| <b>PRODUCT REFERENCE</b> | PTX19021                                |
| <b>CLONALITY</b>         | Monoclonal Antibody                     |
| <b>ISOTYPE</b>           | VHH-8His-Cys-tag                        |
| <b>TARGET</b>            | S Protein                               |

### PRODUCT INFORMATION

---

|                                |   |
|--------------------------------|---|
| <b>HOST SPECIES</b>            | <i>Vicugna pacos (Alpaca)</i>   |
| <b>REACTIVITY</b>              | <i>Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)</i>   |
| <b>FORM</b>                    | Liquid  |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4)  |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**