



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-MERS-CoV S1 N-terminal domain/S1-NTD Monoclonal Antibody, 7D10

OVERVIEW

PRODUCT REFERENCE	PTX19392
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1, S, S1 N-terminal domain, S1-NTD, SARS-like Coronaviruses

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Middle East respiratory syndrome-related coronavirus (isolate UnitedKingdom/H123990006/2012) (MERS-CoV) (Betacoronavirus England 1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.