



Mouse IgG1, kappa (Silence) Isotype Control Antibody (MOPC21)

OVERVIEW

PRODUCT REFERENCE	PTX19426
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	0.01M PBS pH 7.4 and 0.09% sodium azide
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.

