

## PRODUCT DATASHEET

### OVERVIEW

---

|                   |                                       |
|-------------------|---------------------------------------|
| PRODUCT NAME      | Anti-CD107b/LAMP2 Polyclonal Antibody |
| PRODUCT REFERENCE | PTX19577                              |
| CLONALITY         | Polyclonal Antibody                   |
| ISOTYPE           | IgG                                   |

### PRODUCT INFORMATION

---

|                     |                                                                                                                                                                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES      | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                        |
| FORM                | Liquid                                                                                                                                                                                                                                                                                                       |
| BUFFER              |                                                                                                                                                                                                                                                                                                              |
| STABILITY & STORAGE | Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |