



# ProteoGenix

## ProteoGenix S.A.S

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## ProteoGenix INC.

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

|                   |                                            |
|-------------------|--------------------------------------------|
| PRODUCT NAME      | Anti-CD181/CD128/CXCR1 Polyclonal Antibody |
| PRODUCT REFERENCE | PTX19633                                   |
| CLONALITY         | Polyclonal Antibody                        |
| ISOTYPE           | IgG                                        |

### PRODUCT INFORMATION

---

|                     |                                                                                                                                                                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES      | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                        |
| FORM                | Liquid                                                                                                                                                                                                                                                                                                       |
| BUFFER              |                                                                                                                                                                                                                                                                                                              |
| STABILITY & STORAGE | Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY