

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CD130/IL6ST Polyclonal Antibody
PRODUCT REFERENCE	PTX19685
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.