



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-TROP2 & Anti-HSG Bispecific Monoclonal Antibody (TF12)
<b>PRODUCT REFERENCE</b>	PTX25868
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	(Fab) <sub>2</sub> -Fab'

### PRODUCT INFORMATION

---

#### ORIGIN SPECIES

**FORM** Liquid

#### BUFFER

**STABILITY & STORAGE** Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**