



# Proteogenix

## Proteogenix S.A.S

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## Proteogenix INC.

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>PRODUCT NAME</b>      | Anti-Human SERPINE1/PAI-1 (CT140) |
| <b>PRODUCT REFERENCE</b> | PTX26157                          |
| <b>CLONALITY</b>         | Monoclonal Antibody               |
| <b>ISOTYPE</b>           | IgG4, Kappa                       |

### PRODUCT INFORMATION

---

#### ORIGIN SPECIES

**FORM** Liquid

#### BUFFER

**STABILITY & STORAGE** Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**