



Sheep SOCS2WT Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1052
ORIGIN SPECIES	Sheep
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	39,52 kDa
FORM	liquid
BUFFER	Tris 50mM NaCl 150mM, pH9 after refolding
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	W5Q2U4
ALIAS/SYNONYMS	SOCS2WT, suppressor of cytokine signaling 2

