



Trichinella spiralis Adzh68 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1067
ORIGIN SPECIES	Trichinella spiralis
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	44,70 kDa
FORM	liquid
BUFFER	PBS, imidazole 250mM in native conditions. PBS, Urea 4M in denaturing conditions
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	B0F9T9
ALIAS/SYNONYMS	Adzh68-ΔNter, putative serine protease

