



Mouse BRP-39 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1075
ORIGIN SPECIES	House Mouse
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	42,30 kDa
FORM	liquid
BUFFER	PBS, imidazole 10mM, Urea 8M, pH4
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q61362
ALIAS/SYNONYMS	BRP-39, BRP39 protein, Chitinase-3-like protein 1, Chi3l1

