



Rat CFMNter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1079
ORIGIN SPECIES	Rat
FRAGMENT	Partial
TAG	No
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	39,20 kDa
FORM	liquid
BUFFER	TrisHCl 50mM if elution. TrisHCl 50mM, NaCl 150mM, EDTA 1mM, DTT 1mM if cleavage
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6AXS9
ALIAS/SYNONYMS	CFMNter, protein FAM101B, Refilin-B, Regulator of filamin protein B, RefilinB

