



Toxoplasma gondii GRA2(21-185)(RH) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1095
ORIGIN SPECIES	Toxoplasma gondii
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	25,86 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P13404
ALIAS/SYNONYMS	GRA2(21-185)(RH), 28kd antigen

