



Toxoplasma gondii GRA6Nt(41-152)(RH) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1102
ORIGIN SPECIES	Toxoplasma gondii
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	17,46 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM in native conditions. NaH ₂ PO ₄ 100mM, TrisBase 10mM, Urea 8M, imidazole 400mM in denaturig conditions
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q27003
ALIAS/SYNONYMS	GRA6Nt(41-152)(RH), dense granule antigen

