



## Human 3C Protease Recombinant Enzyme (GST tag)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P1106
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Property sequence
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	46,52 kDa
<b>FORM</b>	Frozen
<b>BUFFER</b>	Tris-HC 50mM pH8, NaCl 150mM, EDTA 10mM, DTT 1mM and 20% Glycerol
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	NP_740524.1
<b>ALIAS/SYNONYMS</b>	GST 3C Protease, 3C-Protease

