

HIRA protein, Human HIRA(155-254) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1111
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38,26 kDa
FORM	liquid
BUFFER	TrisHC 50mMl, 10mM glutathione reduced pH8
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P54198
ALIAS/SYNONYMS	HIRA(155-254), HIRA protein