



Human HSP110de9 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1112
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	43,46 kDa
FORM	liquid
BUFFER	PBS, Urea 6M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	/*
ALIAS/SYNONYMS	HSP110de9, heat shock 105kDa protein mutant de9

