

Marinobacter MARHY0478 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1123
ORIGIN SPECIES	Marinobacter
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	46,10 kDa
FORM	liquid
BUFFER	PBS pH 7.4 with imidazole 250mM and Urea 4M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	H8WEC1
ALIAS/SYNONYMS	MARHY0478, long-chain fatty acid transporter, Putative hydrophobic compounds transporter

