



Mycobacterium PptT Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1124
ORIGIN SPECIES	Mycobacterium
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	67,08 kDa
FORM	liquid
BUFFER	Tris 20mM pH7, NaCl 200mM, EDTA 1mM, Maltose 10mM without dialysis. Tris 10mM pH8, NaCl 100mM, glycerol 10% with dialysis
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O33336
ALIAS/SYNONYMS	MBP-PptT, hypothetical protein Rv2794c

