



Ciona intestinalis MYO5 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1127
ORIGIN SPECIES	Ciona intestinalis
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	59,60 kDa
FORM	liquid
BUFFER	PBS pH6-8 containing urea 6M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	63
ALIAS/SYNONYMS	133

