



Trichinella spiralis NBL-1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1128
ORIGIN SPECIES	Trichinella spiralis
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	44,96 kDa
FORM	liquid
BUFFER	NaH ₂ PO ₄ 100mM, TrisHCl 10mM, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9BJL7
ALIAS/SYNONYMS	NBL-1-ΔNter, newborn larvae-specific serine protease SS2

