



Ostreid herpesvirus Putative apoptosis inhibitor ORF87 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1133
ORIGIN SPECIES	Ostreid herpesvirus
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	21,69 kDa
FORM	liquid
BUFFER	PBS, Urea 6M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6R7E2
ALIAS/SYNONYMS	ORF87, ORF87

