

Leishmania P74 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1135
ORIGIN SPECIES	Leishmania
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	50,44 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A4HX73
ALIAS/SYNONYMS	P74, elongation factor 1-alpha

