



Cow PAG1B Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1136
ORIGIN SPECIES	Cow
FRAGMENT	Full-length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	45,67 kDa
FORM	liquid
BUFFER	NaH ₂ PO ₄ 75mM, NaCl 150mM, Urea 6M, pH 6-8.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q29432
ALIAS/SYNONYMS	PAG1B, pregnancy-associated glycoprotein 1 precursor, PAG 1, Pregnancy-specific protein B, PSP-B

