



Plasmodium falciparum PF10Cter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1138
ORIGIN SPECIES	Plasmodium falciparum
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	67,99 kDa
FORM	liquid
BUFFER	TrisHC 50mMl, 10mM reduced glutathion pH8
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8IK14
ALIAS/SYNONYMS	PF10Cter, Tryptophan-rich antigen 3, putative

